

For Research Use Only

LIN7B Polyclonal antibody

Catalog Number: 12732-1-AP



Basic Information

Catalog Number:

12732-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3437

GenBank Accession Number:

BC027618

GeneID (NCBI):

64130

UNIPROT ID:

Q9HAP6

Full Name:

lin-7 homolog B (C. elegans)

Calculated MW:

207 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human heart tissue, A549 cells, human brain tissue, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human heart tissue were subjected to SDS PAGE followed by western blot with 12732-1-AP (LIN7B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.