

For Research Use Only

MKRN1 Polyclonal antibody

Catalog Number: 12773-1-AP



Basic Information

Catalog Number:

12773-1-AP

Size:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3465

GenBank Accession Number:

BC025955

GeneID (NCBI):

23608

UNIPROT ID:

Q9UHC7

Full Name:

makorin ring finger protein 1

Calculated MW:

482 aa, 53 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, Jurkat cells, mouse brain tissue, mouse thymus tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12773-1-AP (MKRN1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.