

For Research Use Only

MKRN1 Polyclonal antibody, PBS Only

Catalog Number:12773-1-PBS



Basic Information

Catalog Number:

12773-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3465

GenBank Accession Number:

BC025955

GeneID (NCBI):

23608

UNIPROT ID:

Q9UHC7

Full Name:

makorin ring finger protein 1

Calculated MW:

482 aa, 53 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12773-1-AP (MKRN1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12773-1-PBS in a different storage buffer formulation.