

For Research Use Only

CNTN4 Polyclonal antibody

Catalog Number: 12777-1-AP



Basic Information

Catalog Number:

12777-1-AP

Size:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3476

GenBank Accession Number:

BC026119

GeneID (NCBI):

152330

UNIPROT ID:

Q8IWW2

Full Name:

contactin 4

Calculated MW:

1026 aa, 113 kDa

Observed MW:

113-130 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12777-1-AP (CNTN4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.