

For Research Use Only

RIC3 Polyclonal antibody

Catalog Number: 12779-1-AP



Basic Information

Catalog Number:

12779-1-AP

Size:

120 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3479

GenBank Accession Number:

BC022455

GeneID (NCBI):

79608

UNIPROT ID:

Q7Z5B4

Full Name:

resistance to inhibitors of
cholinesterase 3 homolog (C.
elegans)

Calculated MW:

368 aa, 41 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

