

For Research Use Only

PAPOLB Polyclonal antibody

Catalog Number: 12821-1-AP



Basic Information

Catalog Number:

12821-1-AP

Size:

220 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3487

GenBank Accession Number:

BC036653

GeneID (NCBI):

56903

UNIPROT ID:

Q9NRJ5

Full Name:

poly(A) polymerase beta (testis specific)

Calculated MW:

71 kDa

Observed MW:

71 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: mouse testis tissue, human kidney tissue, HeLa cells, human brain tissue, human stomach tissue, mouse ovary tissue, human placenta tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 12821-1-AP (PAPOLB antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.