

SECTM1 Polyclonal antibody

Catalog Number: 12865-1-AP

Basic Information

Catalog Number:

12865-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3190

GenBank Accession Number:

BC017716

GeneID (NCBI):

6398

UNIPROT ID:

Q8WVN6

Full Name:

secreted and transmembrane 1

Calculated MW:

248 aa, 27 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

SECTM1, also named as protein K12, is identified by its location just upstream of the CD7 locus. It's a novel molecule with potential importance in hematopoietic and/or immune system processes. The calculated MW of SECTM1 is 27 kDa. It can be cleaved into a 20 kDa soluble form.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

