

For Research Use Only

# SECTM1 Polyclonal antibody

Catalog Number: 12865-1-AP



## Basic Information

<b>Catalog Number:</b> 12865-1-AP	<b>GenBank Accession Number:</b> BC017716	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 133 µg/ml	<b>GeneID (NCBI):</b> 6398	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8WVN6	
<b>Isotype:</b> IgG	<b>Full Name:</b> secreted and transmembrane 1	
<b>Immunogen Catalog Number:</b> AG3190	<b>Calculated MW:</b> 248 aa, 27 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

SECTM1, also named as protein K12, is identified by its location just upstream of the CD7 locus. It's a novel molecule with potential importance in hematopoietic and/or immune system processes. The calculated MW of SECTM1 is 27 kDa. It can be cleaved into a 20 kDa soluble form.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data