

For Research Use Only

# RPP38 Polyclonal antibody, PBS Only

Catalog Number:12874-1-PBS



## Basic Information

Catalog Number:

12874-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3552

GenBank Accession Number:

BC029494

GeneID (NCBI):

10557

UNIPROT ID:

P78345

Full Name:

ribonuclease P/MRP 38kDa subunit

Calculated MW:

38 kDa

Observed MW:

38 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



rat kidney tissue were subjected to SDS PAGE followed by western blot with 12874-1-AP (RPP38 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12874-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using RPP38 antibody (12874-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 12874-1-PBS in a different storage buffer formulation.