

For Research Use Only

LILRB5 Polyclonal antibody

Catalog Number: 12881-1-AP



Basic Information

Catalog Number:

12881-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3886

GenBank Accession Number:

BC025704

GeneID (NCBI):

10990

UNIPROT ID:

O75023

Full Name:

leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5

Calculated MW:

590 aa, 64 kDa

Observed MW:

64 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : COLO 320 cells, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12881-1-AP (ILRB5 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.