

For Research Use Only

PDE4D Polyclonal antibody

Catalog Number: 12918-1-AP

Featured Product

24 Publications



Basic Information

Catalog Number:

12918-1-AP

Concentration:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3588

GenBank Accession Number:

BC036319

GeneID (NCBI):

5144

UNIPROT ID:

Q08499

Full Name:

phosphodiesterase 4D, cAMP-specific
(phosphodiesterase E3 dunce
homolog, Drosophila)

Calculated MW:

809 aa, 91 kDa

Observed MW:

91 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IHC 1:50-1:500

IF/ICC 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat, rabbit

Positive Controls:

WB : K-562 cells, mouse brain tissue, mouse pancreas
tissue

IHC : mouse colon tissue,

IF/ICC : HeLa cells,

**Note-IHC: suggested antigen retrieval with
TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Background Information

PDE4D(cAMP-specific 3',5'-cyclic phosphodiesterase 4D) is also named as DPDE3, PDE43 and belongs to the cyclic nucleotide phosphodiesterase family and PDE4 subfamily. PDE4D subtype may be of particular importance as an antidepressant target in that it is regulated by repeated treatment with both norepinephrine and serotonin reuptake inhibitors as well as by the PDE4 inhibitor rolipram(PMID:16730340). It has some isoforms and all the PDE4D isoforms are expressed in healthy adult human hippocampus and in hippocampus from a patient with advanced Alzheimer disease(PMID:17353396).

Notable Publications

Author	Pubmed ID	Journal	Application
Keita Tsujimura	26344767	Cell Rep	WB
Fude Liu	31772658	J Cancer	WB
Xiang Tao	34826603	Pharmacol Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




K562 cell lysate were subjected to SDS PAGE followed by western blot with 12918-1-AP (PDE4D antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using 12918-1-AP (PDE4D Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using 12918-1-AP (PDE4D Antibody) at dilution of 1:100 (under 40x lens).



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using PDE4D antibody (12918-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L).