

For Research Use Only

GCAP1 Polyclonal antibody

Catalog Number: 12944-1-AP



Basic Information

Catalog Number:

12944-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4017

GenBank Accession Number:

BC031663

GeneID (NCBI):

2978

UNIPROT ID:

P43080

Full Name:

guanylate cyclase activator 1A
(retina)

Calculated MW:

23 kDa

Observed MW:

23-24 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Y79 cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Y79 cells were subjected to SDS PAGE followed by western blot with 12944-1-AP (GCAP1 antibody) at dilution of 1:200 incubated at room temperature for 1.5 hours.