

For Research Use Only

# IMPDH2 Polyclonal antibody, PBS Only

Catalog Number: 12948-1-PBS

Featured Product



## Basic Information

Catalog Number:

12948-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3632

GenBank Accession Number:

BC012840

GeneID (NCBI):

3615

UNIPROT ID:

P12268

Full Name:

IMP (inosine monophosphate)  
dehydrogenase 2

Calculated MW:

514 aa, 56 kDa

Observed MW:

56 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, IF/ICC, IP, ELISA

Species Specificity:

human, mouse, rat

## Background Information

IMPDH2, also named as IMPD2, belongs to the IMPDH/GMPR family. It is a rate limiting enzyme in the de novo synthesis of guanine nucleotides and therefore is involved in the regulation of cell growth. It may also have a role in the development of malignancy and the growth progression of some tumors. IMPDH2 is a key enzyme in the purine nucleotide biosynthetic pathway and constitutes a pivotal biological target for immunosuppressant and antiviral drugs. (PMID:21181270) IMPDH2 is highly enriched in RR(filamentous rods and rings). This antibody may have cross-reaction with IMPDH1 due to the high homology.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 12948-1-AP (IMPDH2 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 12948-1-PBS in a different storage buffer formulation.




WB result of IMPDH2 antibody (12948-1-AP; 1:12000; incubated at room temperature for 1.5 hours) with sh-Control and sh-IMPDH2 transfected A549 cells. This data was developed using the same antibody clone with 12948-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 12948-1-AP (IMPDH2 antibody) at dilution of 1:18000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 12948-1-PBS in a different storage buffer formulation.



IP result of anti-IMPDH2 (IP:12948-1-AP, 4ug; Detection:12948-1-AP 1:20000) with HeLa cells lysate 1320 ug. This data was developed using the same antibody clone with 12948-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using IMPDH2 antibody (12948-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-phalloidin (red). This data was developed using the same antibody clone with 12948-1-PBS in a different storage buffer formulation.