

For Research Use Only

# JAM2 Polyclonal antibody

Catalog Number: 12972-1-AP

2 Publications



## Basic Information

Catalog Number:

12972-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3431

GenBank Accession Number:

BC017779

GeneID (NCBI):

58494

UNIPROT ID:

P57087

Full Name:

junctional adhesion molecule 2

Calculated MW:

298 aa, 33 kDa

Observed MW:

33 kDa, 43 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : human placenta tissue, human heart tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Hiroki Takeuchi	34448537	Cell Microbiol	WB
Guang Xu	28081450	Diabetes Res Clin Pract	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 12972-1-AP (JAM2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.