

For Research Use Only

JAM2 Polyclonal antibody

Catalog Number: 12972-1-AP **2 Publications**



Basic Information

Catalog Number: 12972-1-AP	GenBank Accession Number: BC017779	Purification Method: Antigen affinity purification
Size: 600 µg/ml	GeneID (NCBI): 58494	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: P57087	
Isotype: IgG	Full Name: junctional adhesion molecule 2	
Immunogen Catalog Number: AG3431	Calculated MW: 298 aa, 33 kDa	
	Observed MW: 33 kDa, 43 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human placenta tissue, human heart tissue
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hiroki Takeuchi	34448537	Cell Microbiol	WB
Guang Xu	28081450	Diabetes Res Clin Pract	WB


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 12972-1-AP (JAM2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.