

For Research Use Only

ARMER Polyclonal antibody

Catalog Number: 12988-1-AP



Basic Information

Catalog Number:

12988-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3639

GenBank Accession Number:

BC010281

GeneID (NCBI):

23204

UNIPROT ID:

Q15041

Full Name:

ADP-ribosylation factor-like 6
interacting protein 1

Calculated MW:

203 aa, 23 kDa

Observed MW:

22-24 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse lung tissue, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 12988-1-AP (ARMER antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.