

For Research Use Only

SIVA1 Polyclonal antibody

Catalog Number: 12992-1-AP



Basic Information

Catalog Number:

12992-1-AP

Size:

330 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3654

GenBank Accession Number:

BC034562

GeneID (NCBI):

10572

UNIPROT ID:

O15304

Full Name:

SIVA1, apoptosis-inducing factor

Calculated MW:

175 aa, 19 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

Background Information

SIVA1, also named as SIVA, CD27BP, is an apoptosis regulatory protein. It induces the CD27-mediated apoptosis. SIVA is required for its binding to BCL-XL and sensitizing cells to UV radiation.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human ovary tumor using 12992-1-AP (SIVA1 antibody) at dilution of 1:100 (under 10x lens).