

For Research Use Only

p18 INK4c Polyclonal antibody

Catalog Number: 12993-1-AP



Basic Information

Catalog Number:

12993-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3655

GenBank Accession Number:

BC016173

GeneID (NCBI):

1031

UNIPROT ID:

P42773

Full Name:

cyclin-dependent kinase inhibitor 2C
(p18, inhibits CDK4)

Calculated MW:

168 aa, 18 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer using 12993-1-AP (p18 INK4c antibody) at dilution of 1:100 (under 10x lens).