

For Research Use Only

CEBPG Polyclonal antibody

Catalog Number: 12997-1-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

12997-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3663

GenBank Accession Number:

BC013128

GeneID (NCBI):

1054

UNIPROT ID:

P53567

Full Name:

CCAAT/enhancer binding protein (C/EBP), gamma

Calculated MW:

150 aa, 16 kDa

Observed MW:

16-18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:250-1:1000

IF/ICC 1:1000-1:4000

Applications

Tested Applications:

IHC, IF/ICC, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : mouse liver tissue, mouse colon tissue, rat liver tissue

IF/ICC : MCF-7 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yongsheng Huang	33163273	Am J Cancer Res	WB,IHC
Feng Tian	35767995	Neuron	IHC
Jiahong Tan	39238655	Cancer Inform	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 12997-1-AP (CEBPG antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 12997-1-AP (CEBPG antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CEBPG antibody (12997-1-AP) at dilution of 1:2000 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using 12997-1-AP (CEBPG antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using 12997-1-AP (CEBPG antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).