

For Research Use Only

# ZMYND11 Polyclonal antibody

Catalog Number: **13012-1-AP**



## Basic Information

Catalog Number:

**13012-1-AP**

Size:

**N/A**

Source:

**Rabbit**

Isotype:

**IgG**

Immunogen Catalog Number:

**AG4016**

GenBank Accession Number:

**BC034784**

GeneID (NCBI):

**10771**

UNIPROT ID:

**Q15326**

Full Name:

**zinc finger, MYND domain containing**

**11**

Calculated MW:

**527 aa, 62 kDa**

Purification Method:

**Antigen affinity purification**

## Applications

Tested Applications:

**ELISA**

Species Specificity:

**human, mouse, rat**

## Background Information

### Storage

Storage:

**Store at -20°C. Stable for one year after shipment.**

Storage Buffer:

**PBS with 0.02% sodium azide and 50% glycerol pH 7.3.**

**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

