

For Research Use Only

# PLEKHA5 Polyclonal antibody

Catalog Number: 13045-1-AP



## Basic Information

Catalog Number:

13045-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3575

GenBank Accession Number:

BC013133

GeneID (NCBI):

54477

UNIPROT ID:

Q9HAU0

Full Name:

pleckstrin homology domain  
containing, family A member 5

Calculated MW:

1116 aa, 127 kDa

Observed MW:

127 kDa, 74 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: mouse testis tissue, human heart tissue, mouse  
kidney tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13045-1-AP (PLEKHA5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.