

For Research Use Only

PLEKHA5 Polyclonal antibody

Catalog Number: 13045-1-AP



Basic Information

Catalog Number:

13045-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3575

GenBank Accession Number:

BC013133

GeneID (NCBI):

54477

UNIPROT ID:

Q9HAU0

Full Name:

pleckstrin homology domain
containing, family A member 5

Calculated MW:

1116 aa, 127 kDa

Observed MW:

127 kDa, 74 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, human heart tissue, mouse
kidney tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13045-1-AP (PLEKHA5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.