

For Research Use Only

# MCFD2 Polyclonal antibody

Catalog Number: 13068-1-AP



## Basic Information

Catalog Number:

13068-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3768

GenBank Accession Number:

BC040357

GeneID (NCBI):

90411

UNIPROT ID:

Q8NI22

Full Name:

multiple coagulation factor  
deficiency 2

Calculated MW:

146 aa, 16 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF 1:20-1:200

## Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

IF : HepG2 cells, HeLa cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of HepG2 cells using 13068-1-AP (MCFD2 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).