

For Research Use Only

# NUP93 Polyclonal antibody

Catalog Number: 13077-1-AP

Featured Product

6 Publications



## Basic Information

Catalog Number:

13077-1-AP

Concentration:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3735

GenBank Accession Number:

BC034346

GeneID (NCBI):

9688

UNIPROT ID:

Q8N1F7

Full Name:

nucleoporin 93kDa

Calculated MW:

819 aa, 93 kDa

Observed MW:

85-93 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, A431 cells, HEK-293 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhixiong Liu	30876848	Neuron	WB
Yuri Furusawa	30087672	Front Microbiol	WB
Bin Chen	35152263	Oncogene	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13077-1-AP (NUP93 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.