

For Research Use Only

RHOXF2 Polyclonal antibody

Catalog Number: 13088-1-AP **1 Publications**



Basic Information

Catalog Number: 13088-1-AP	GenBank Accession Number: BC021719	Purification Method: Antigen affinity purification
Size: 1000 µg/ml	GeneID (NCBI): 84528	Recommended Dilutions: WB 1:500-1:2400
Source: Rabbit	UNIPROT ID: Q9BQY4	
Isotype: IgG	Full Name: Rhox homeobox family, member 2	
Immunogen Catalog Number: AG3602	Calculated MW: 32 kDa	
	Observed MW: 35-38 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : K-562 cells,
Cited Applications: WB, IF, IHC	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jennifer Borgmann	28171660	Hum Mol Genet	WB, IHC, IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 13088-1-AP (RHOXF2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.