

For Research Use Only

FGF13 Polyclonal antibody

Catalog Number: 13201-1-AP **3 Publications**



Basic Information

Catalog Number: 13201-1-AP	GenBank Accession Number: BC034340	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 2258	
Source: Rabbit	UNIPROT ID: Q92913	
Isotype: IgG	Full Name: fibroblast growth factor 13	
Immunogen Catalog Number: AG3942	Calculated MW: 28 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB, IF

Species Specificity:
human, mouse, rat

Cited Species:
human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhang Xiujun X	23554914	PLoS One	WB
Xiaojing Yue	29405966	Placenta	WB,IF
J-J Song	30720165	Eur Rev Med Pharmacol Sci	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data