

For Research Use Only

# FGF13 Polyclonal antibody

Catalog Number: 13201-1-AP

3 Publications



## Basic Information

Catalog Number:

13201-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3942

GenBank Accession Number:

BC034340

GeneID (NCBI):

2258

UNIPROT ID:

Q92913

Full Name:

fibroblast growth factor 13

Calculated MW:

28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhang Xiujun X	23554914	PLoS One	WB
Xiaojing Yue	29405966	Placenta	WB,IF
J-J Song	30720165	Eur Rev Med Pharmacol Sci	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

