

For Research Use Only

Granzyme M Polyclonal antibody

Catalog Number: 13208-1-AP



Basic Information

Catalog Number:

13208-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3982

GenBank Accession Number:

BC025701

GeneID (NCBI):

3004

UNIPROT ID:

P51124

Full Name:

granzyme M (Lymphocyte met-ase 1)

Calculated MW:

257 aa, 27 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse lymph tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lymph tissue were subjected to SDS PAGE followed by western blot with 13208-1-AP (GZMM antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.