

For Research Use Only

# ATCAY Polyclonal antibody

Catalog Number: 13257-2-AP



## Basic Information

Catalog Number:

13257-2-AP

Concentration:

280 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3658

GenBank Accession Number:

BC026217

GeneID (NCBI):

85300

UNIPROT ID:

Q86WG3

Full Name:

ataxia, cerebellar, Cayman type

Calculated MW:

371 aa, 42 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13257-2-AP (ATCAY antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.