

For Research Use Only

PLZF Polyclonal antibody

Catalog Number: 13259-1-AP

1 Publications



Basic Information

Catalog Number:

13259-1-AP

Size:

520 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3783

GenBank Accession Number:

BC026902

GeneID (NCBI):

7704

UNIPROT ID:

Q05516

Full Name:

zinc finger and BTB domain containing 16

Calculated MW:

673 aa, 74 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

goat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Huanshan He	38167777	J Anim Sci	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

