

For Research Use Only

EDN2 Polyclonal antibody

Catalog Number: 13272-1-AP



Basic Information

Catalog Number:

13272-1-AP

Size:

310 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3969

GenBank Accession Number:

BC034393

GeneID (NCBI):

1907

UNIPROT ID:

P20800

Full Name:

endothelin 2

Calculated MW:

178 aa, 20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung cancer using 13272-1-AP (EDN2 antibody) at dilution of 1:100 (under 40x lens).