

For Research Use Only

UNC5CL Polyclonal antibody

Catalog Number: 13322-1-AP



Basic Information

Catalog Number: 13322-1-AP	GenBank Accession Number: BC035284	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 222643	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q8IV45	
Isotype: IgG	Full Name: unc-5 homolog C (C. elegans)-like	
Immunogen Catalog Number: AG3970	Calculated MW: 518 aa, 58 kDa	
	Observed MW: 69 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse colon tissue,
Species Specificity: human, mouse	

Background Information

Unc5CL is a factor involved in epithelial inflammation and immunity.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse colon tissue were subjected to SDS PAGE followed by western blot with 13322-1-AP (UNC5CL Antibody) at dilution of 1:300 incubated at 4 degree celsius over night.