

For Research Use Only

PI3K p101 Polyclonal antibody

Catalog Number: 13329-1-AP **8 Publications**



Basic Information

Catalog Number:

13329-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4137

GenBank Accession Number:

BC028212

GeneID (NCBI):

23533

UNIPROT ID:

Q8WYR1

Full Name:

phosphoinositide-3-kinase,
regulatory subunit 5

Calculated MW:

880 aa, 97 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rabbit

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yiyi Jin	29285087	Exp Ther Med	WB
Xiaoqin Zhu	34887935	Evid Based Complement Alternat Med	WB
Hong-Yan Lin	32175001	Chin Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data