

For Research Use Only

# LENG8 Polyclonal antibody, PBS Only

Catalog Number:13336-1-PBS



## Basic Information

<b>Catalog Number:</b> 13336-1-PBS	<b>GenBank Accession Number:</b> BC028048	<b>Purification Method:</b> Antigen affinity purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 114823	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> Q96PV6	
<b>Immunogen Catalog Number:</b> AG4153	<b>Full Name:</b> Leukocyte receptor cluster (LRC) member 8	
	<b>Calculated MW:</b> 800 aa, 88 kDa	
	<b>Observed MW:</b> 86-88 kDa	

## Applications

**Tested Applications:**  
WB, FC (Intra), Indirect ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

LENG8 is the core factor of the key RNA quality control system discovered in recent years, which plays a goalkeeper role in the regulation of eukaryotic gene expression. LENG8 promotes the recruitment of PAXT complex and RNA exonuclease, and mediates the degradation of abnormal RNA. It is a multifunctional evolutionarily conserved protein, which plays a key role from yeast to human beings. It is an important node molecule connecting gene post-transcriptional processing, RNA quality control and cell energy metabolism.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



A2780 cells were subjected to SDS PAGE followed by western blot with 13336-1-AP (LENG8 antibody) at dilution of 1:1200 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 13336-1-PBS in a different storage buffer formulation.



Jurkat cells were subjected to SDS PAGE followed by western blot with 13336-1-AP (LENG8 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 13336-1-PBS in a different storage buffer formulation.



1x10<sup>6</sup> HEK-293 cells were intracellularly stained with 0.25 ug LENG8 Polyclonal antibody (13336-1-AP) and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011). This data was developed using the same antibody clone with 13336-1-PBS in a different storage buffer formulation.