

For Research Use Only

RPP14 Polyclonal antibody

Catalog Number: 13339-1-AP



Basic Information

Catalog Number:

13339-1-AP

Size:

190 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4159

GenBank Accession Number:

BC012017

GeneID (NCBI):

11102

UNIPROT ID:

O95059

Full Name:

ribonuclease P/MRP 14kDa subunit

Calculated MW:

124 aa, 14 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

RPP14 is a subunit of human ribonuclease P. RPP14 may be part of an exosome activity. This antibody is specific to RPP14.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

