

For Research Use Only

KCNRG Polyclonal antibody

Catalog Number: 13344-1-AP



Basic Information

Catalog Number:

13344-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4176

GenBank Accession Number:

BC032343

GeneID (NCBI):

283518

UNIPROT ID:

Q8N5I3

Full Name:

potassium channel regulator

Calculated MW:

272 aa, 31 kDa

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse lung tissue, rat heart tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13344-1-AP (KCNRG antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.