

For Research Use Only

# NLRC4 Polyclonal antibody

Catalog Number: 13390-1-AP

1 Publications



## Basic Information

Catalog Number:

13390-1-AP

Size:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4048

GenBank Accession Number:

BC031555

GeneID (NCBI):

58484

UNIPROT ID:

Q9NPP4

Full Name:

NLR family, CARD domain containing

4

Calculated MW:

1024 aa, 116 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Shuai Peng	39060654	Tissue Engineering And Regenerative Medicine	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

