

For Research Use Only

PAQR4 Polyclonal antibody

Catalog Number: 13401-1-AP

6 Publications



Basic Information

Catalog Number:

13401-1-AP

Size:

153 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4241

GenBank Accession Number:

BC033703

GeneID (NCBI):

124222

UNIPROT ID:

Q8N4S7

Full Name:

progesterone and adiponectin receptor family member IV

Calculated MW:

273 aa, 29 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : HeLa cells, mouse liver tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiaobao Shi	36241701	Cancer Gene Ther	WB
Gang Zhao	34718369	Acta Biochim Biophys Sin (Shanghai)	WB
Wenbin Chen	33490086	Front Cell Dev Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 13401-1-AP (PAQR4 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.