

For Research Use Only

ACYP1 Polyclonal antibody

Catalog Number: 13433-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

13433-1-AP

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4237

GenBank Accession Number:

BC035568

GeneID (NCBI):

97

UNIPROT ID:

P07311

Full Name:

acylphosphatase 1, erythrocyte (common) type

Calculated MW:

99 aa, 11 kDa

Observed MW:

11 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB: mouse testis tissue, rat testis tissue, human brain tissue

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|------------|-----------|-------------------|-------------|
| Shuai Wang | 37210811 | Drug Resist Updat | WB,IHC,IP |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13433-1-AP (ACYP1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.