

For Research Use Only

RASL11B Polyclonal antibody, PBS Only

Catalog Number: 13476-1-PBS



Basic Information

Catalog Number:

13476-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4277

GenBank Accession Number:

BC025694

GeneID (NCBI):

65997

UNIPROT ID:

Q9BPW5

Full Name:

RAS-like, family 11, member B

Calculated MW:

248 aa, 18 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data