

For Research Use Only

RHEBL1 Polyclonal antibody

Catalog Number: 13478-1-AP



Basic Information

Catalog Number:

13478-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4286

GenBank Accession Number:

BC027482

GeneID (NCBI):

121268

UNIPROT ID:

Q8TAI7

Full Name:

Ras homolog enriched in brain like 1

Calculated MW:

183 aa, 21 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF-P, ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

