

For Research Use Only

CDH23 Polyclonal antibody

Catalog Number: 13496-1-AP

1 Publications



Basic Information

Catalog Number:

13496-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4322

GenBank Accession Number:

BC032581

GeneID (NCBI):

64072

UNIPROT ID:

Q9H251

Full Name:

cadherin-like 23

Calculated MW:

3354 aa, 369 kDa

Observed MW:

59 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : A431 cells, Y79 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Fabian Jansen	26620972	Hum Mol Genet	WB, IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 13496-1-AP (CDH23 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.