

For Research Use Only

cyclin G2 Polyclonal antibody

Catalog Number: 13535-1-AP

1 Publications



Basic Information

Catalog Number:

13535-1-AP

Size:

100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4449

GenBank Accession Number:

BC032518

GeneID (NCBI):

901

UNIPROT ID:

Q16589

Full Name:

cyclin G2

Calculated MW:

344 aa, 39 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, WB, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jinlan Gao	30547803	J Exp Clin Cancer Res	WB, IHC, IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

