

For Research Use Only

JAK3 Polyclonal antibody

Catalog Number: 13551-1-AP **1 Publications**



Basic Information

Catalog Number: 13551-1-AP	GenBank Accession Number: BC028068	Purification Method: Antigen affinity purification
Size: 70 µg/ml	GeneID (NCBI): 3718	
Source: Rabbit	UNIPROT ID: P52333	
Isotype: IgG	Full Name: Janus kinase 3 (a protein tyrosine kinase, Leukocyte)	
Immunogen Catalog Number: AG4298	Calculated MW: 1124 aa, 125 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
Human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qiqiang Zhang	33228635	BMC Complement Med Ther	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data