

For Research Use Only

# RAD51L1 Polyclonal antibody

Catalog Number: 13555-1-AP

1 Publications



## Basic Information

Catalog Number:

13555-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4412

GenBank Accession Number:

BC030219

GeneID (NCBI):

5890

UNIPROT ID:

O15315

Full Name:

RAD51-like 1 (S. cerevisiae)

Calculated MW:

384 aa, 42 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Liming Huang	37858068	Cancer Cell Int	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 13555-1-AP (RAD51L1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.