

For Research Use Only

HCFC2 Polyclonal antibody

Catalog Number: 13567-1-AP



Basic Information

Catalog Number:

13567-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4461

GenBank Accession Number:

BC033799

GeneID (NCBI):

29915

UNIPROT ID:

Q9Y5Z7

Full Name:

host cell factor C2

Calculated MW:

792 aa, 87 kDa

Observed MW:

87 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13567-1-AP (HCFC2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.