

For Research Use Only

RASSF4 Polyclonal antibody

Catalog Number: **13575-1-AP**



Basic Information

Catalog Number:

13575-1-AP

Size:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4475

GenBank Accession Number:

BC032593

GeneID (NCBI):

83937

UNIPROT ID:

Q9H2L5

Full Name:

**Ras association (RalGDS/AF-6)
domain family member 4**

Calculated MW:

321 aa, 37 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

