

For Research Use Only

ST3GAL2 Polyclonal antibody

Catalog Number: 13620-1-AP



Basic Information

Catalog Number: 13620-1-AP	GenBank Accession Number: BC036777	Purification Method: Antigen affinity purification
Size: 240 µg/ml	GeneID (NCBI): 6483	Recommended Dilutions: WB 1:500-1:1000 IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q16842	
Isotype: IgG	Full Name: ST3 beta-galactoside alpha-2,3-sialyltransferase 2	
Immunogen Catalog Number: AG4528	Calculated MW: 350 aa, 40 kDa	
	Observed MW: 55 kDa	

Applications

Tested Applications:

IHC, WB, ELISA

Species Specificity:

human, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: human colon tissue, HL-60 cells, human heart tissue, human liver tissue, human spleen tissue, K-562 cells

IHC: rat heart tissue, human skin tissue

Background Information

Beta-galactoside alpha-2,3-sialyltransferase 2 (ST3GAL2) is a sialyltransferase (ST) that transfers sialic acid from the activated donor CMP-sialic acid to the glycan terminus of a glycolipid or glycoprotein. ST3GAL2 contains two putative N-glycosylation sites (Asn92 and Asn211), which are mainly responsible for GD1a and GT1b ganglioside biosynthesis in the brain. ST3GAL2 belongs to the sialyltransferase family, which mediates the terminal modification of glycoproteins and glycolipids. ST3GAL2 has been found to play a role in obesity, aging, and malignant diseases (PMID: 34319453). The calculated molecular weight of ST3GAL2 is 40 kDa, but the observed is 55 kDa, which may be the glycosylation form of ST3GAL2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human colon tissue were subjected to SDS PAGE followed by western blot with 13620-1-AP (ST3GAL2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded rat heart tissue slide using 13620-1-AP (ST3GAL2 antibody) at dilution of 1:100 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).