

For Research Use Only

TPPP2 Polyclonal antibody

Catalog Number: 13633-1-AP **1 Publications**



Basic Information

Catalog Number:

13633-1-AP

Size:

1800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4561

GenBank Accession Number:

BC038970

GeneID (NCBI):

122664

UNIPROT ID:

P59282

Full Name:

tubulin polymerization-promoting
protein family member 2

Calculated MW:

170 aa, 18 kDa

Observed MW:

18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yiling Zhou	39704931	Sci China Life Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13633-1-AP (TPPP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.