

For Research Use Only

# MAK Polyclonal antibody

Catalog Number: 13638-1-AP **1 Publications**



## Basic Information

Catalog Number:

13638-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4571

GenBank Accession Number:

BC039825

GeneID (NCBI):

4117

UNIPROT ID:

P20794

Full Name:

male germ cell-associated kinase

Calculated MW:

623 aa, 71 kDa

Observed MW:

60-62 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Positive Controls:

WB : PC-3 cells, HeLa cells

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Ran Li	28073914	J Biol Chem	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 13638-1-AP (MAK antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.